

Validation Kit

For MicroBio Bioaerosol Samplers



Developed specifically for use with the 100 litre per minute MicroBio MB1 and MB2 MicroBio bioaerosol samplers, this validation kit could not be simpler to use.

This Validation Kit will ensure your MicroBio Air Sampler is performing perfectly and satisfies industry regulators and standards.

Regular validation is required for regulatory compliance during normal use in addition to sending the device away for routine calibration. This is not a substitute for annual calibration of the MicroBio Air Samplers. This service can be undertaken through your local distributor or sending your device directly to Cantium Scientific Limited.

The kit requires calibration only every 5 years and is manufactured to exacting dimensions, referenced to national standards. This will keep your operating costs to a minimum.



CANTIUM SCIENTIFIC
ENGINEERING FOR SCIENCE

Key Features

- ▶ Calibration traceable to National Standards.
- ▶ Lightweight and easy to use.
- ▶ Does not require batteries or power of any kind.
- ▶ Each unit referenced against standard prior to shipment.
- ▶ Total weight 1.6kg (including case and all accessories).
- ▶ Simple visual indication of the calibration "hot-spot" means minimal operator training is required.
- ▶ Designed and manufactured in the UK under ISO 9001:2008 quality systems.

Supplied with

- ▶ Rugged foam lined plastic case.
- ▶ Flow rig, 90 mm Petri dish, 55 mm contact plate and miniature screwdriver.
- ▶ Calibration Certificate.

Dimensions

- ▶ Case 350 mm x 300 mm x 90 mm